



DATA SHEET

Dil-HDL

High Density Lipoprotein labeled with 1,1'-dioctadecyl-3,3,3',3'-tetramethyl-indocarbocyanine perchlorate

Catalog No: YB-0017

Concentration: >2.0mg/ml (Corrected by Cholesterol in Lipoprotein, HDL-C:4.88mmol/L, Protein assayed by BCA:7.60mg/ml)

Quantity: >2.0mg(1ml)/vial

Lot No: lot specific

Absorbance Ratio: $\frac{\text{Dil}}{\text{Protein}} = \frac{555\text{nm}}{275\text{nm}} = 1.82$

Product Preparation:

Purified HDL is labeled with the fluorescent probe, Dil, and reisolated by ultracentrifugation (1.063 ~1.21g/mL). The resultant product is exhaustively dialyzed against phosphate buffered saline, sterilized by membrane filtration and then aseptically packaged under nitrogen in a solution containing phosphate-buffered saline at pH 7.4 and 0.2 mM EDTA-Na₂.

Storage & Stability:

Dil-HDL should be kept sterile at 2-8°C. **NEVER FREEZE**. The stability of this product is in the 6 week range after receipt. Clarify by centrifugation if needed.

Typical Lipoprotein Labeling Protocol

1. Aseptically dilute the Dil-HDL to 10-40µg/ml in your culture media.
2. Add to live cells and incubate for 4-6 hours at 37°C.
3. Remove media containing Dil-HDL from your culture.
4. Wash cells several times with probe-free media.
5. A. Fluorescence Microscopy:

Visualize using standard rhodamine excitation: emission filters (or suggested wavelengths excitation:emission at 554nm:571nm or near). If fixation is desired use 3% formaldehyde in PBS. (Never use methanol or acetone fixation - Dil is soluble in organic solvents). Note: A positive culture must be stained for comparison purposes.

B. Cell Sorting:

Label as in steps 1-5. Trypsinize or treat cultures with EDTA to produce a single cell suspension. Use labeled pure cultures of positive and negative cell types to set gates on the cell sorter.

Suggested Wavelengths for Cell Sorting: Excitation: 554nm
Emission: 571nm

Fixation and Mounting of Dil Labeled Cells

1. Wash 3 times in PBS.
2. Fix in 4% formaldehyde/PBS for 20 minutes at room temperature.
3. Rinse 5 seconds in distilled water at room temperature.
4. Drain liquid onto chem-wipe.
5. Invert cover, slip on a drop of 90% Glycerol and 10% PBS onto a microscope slide.
6. Seal with Kroenigs wax, also known as cover glass cement (Pfaltz & Bauer, Waterbury, CT 06708). Do not use nail polish. Store at -20°C.

***Special Note:** HDL products have a natural tendency to aggregate. Aggregates of this product can interfere with its use. To clarify these aggregates out, simply spin in a microfuge for 2 minutes.

FOR RESEARCH USE ONLY